



## PROTOCOL: Direct PCR from Blood Dried onto Filter Paper

v. 20210810

### Sample:

- Blood dried onto Whatman Filter paper
- Other samples may be suitable (must be verified by user)

### Materials:

- Entopsis:
  - PCRopsis™ Elution Buffer
  - PCRopsis™ Reagent RVD with RVD Enhancer
- PCR master mix
- Primers
- Test sample
- Hole puncher
- Thin-walled PCR tubes
- qPCR Thermocycler

### Methods:

1. Create hole punches of your test sample using a standard, sterile hole puncher
2. Place 1 hole punched filter paper sample into a 1.5 mL tube
3. Add ~100  $\mu$ L of PCRopsis™ Elution Buffer to the tube
4. Vortex thoroughly for 1~3 minutes to elute material off the filter paper
5. Mix 20  $\mu$ L PCRopsis™ Reagent RVD with RVD Enhancer with 20  $\mu$ L of eluted sample in a thin-walled tube (0.2 ~ 0.6 mL) and cap tubes
6. Heat mixture of reagent + blood sample for 10 minutes at 95°C and let cool at room temperature for ~10 seconds before continuing
7. Mix thoroughly
8. Add 5  $\mu$ L of heated sample (top, clear phase of heated sample) to 15  $\mu$ L of qPCR mixture.
  - a. qPCR mixture:
    - i. 5x PCR Master Mix: 4  $\mu$ L
    - ii. Forward primer (10 pmol/  $\mu$ L): 1  $\mu$ L
    - iii. Reverse primer (10 pmol/  $\mu$ L): 1  $\mu$ L
    - iv. Probe (10 pmol/  $\mu$ L): 0.5  $\mu$ L
    - v. Nuclease-Free Water: 8.5  $\mu$ L
9. Run samples on qPCR Thermocycler for 45 cycles.
  - a. DNA Amplification:
    - i. 95°C 5 minutes (initial denaturation)
      1. 95°C 30 seconds
      2. 55°C 30 seconds
      3. 72°C 30 seconds
    - ii. 72°C 60 seconds (final extension)
    - iii. 4°C hold