



PROTOCOL: Direct PCR from Urine Samples

v. 20210810

Sample:

- Urine sample
 - bacterial cells were pelleted at $>1,400\times g$ for ~10 minutes
 - most of the urine supernatant was removed following centrifugation, leaving only 100 ~ 500 μL of urine over the pellet

Materials:

- Entopsis:
 - PCR_{opsis}[™] Support
 - PCR_{opsis}[™] Reagent RVD with RVD Enhancer
- PCR master mix
- Primers
- Test sample
- Thin-walled PCR tubes
- qPCR Thermocycler

Methods:

1. Add 2 μL of PCR_{opsis}[™] Support to 1 mL PCR_{opsis}[™] Reagent RVD with RVD Enhancer, and mix thoroughly
 - a. This mixture is stable for 24 hours
2. Thoroughly mix 20 μL of your sample + 20 μL PCR_{opsis}[™] Reagent RVD with RVD Enhancer with PCR_{opsis}[™] Support in a thin-walled tube (0.2 ~ 0.6 mL) and cap tube
3. Heat mixture of reagent + urine sample for 10 minutes at 95°C and let cool at room temperature for ~10 seconds before continuing
4. Mix thoroughly
5. Add 5 μL of heated sample to 15 μL of qPCR mixture.
 - a. qPCR mixture:
 - i. 5x PCR Master Mix: 4 μL
 - ii. Forward primer (10 pmol/ μL): 1 μL
 - iii. Reverse primer (10 pmol/ μL): 1 μL
 - iv. Probe (10 pmol/ μL): 0.5 μL
 - v. Nuclease-Free Water: 8.5 μL
6. Run samples on qPCR Thermocycler for 45 cycles.
 - a. DNA Amplification:
 - i. 95°C 5 minutes (initial denaturation)
 1. 95°C 30 seconds
 2. 55°C 30 seconds
 3. 72°C 30 seconds
 - ii. 72°C 60 seconds (final extension)
 - iii. 4°C hold