

# PCRopsis™ Reagent RVD-E

## For In Vitro Diagnostic Use

REF #: 7833025, 7833100 & 7833500

REF #: 7833001 (for validation purposes only)

Store at room temperature

### **INTENDED USE**

PCRopsis™ Reagent RVD-E is intended for extraction-free amplification of RNA and DNA from specimens on swabs, without the need for transport mediums.

## 01 INTRODUCTION

PCR<sub>opsis</sub><sup>™</sup> Reagent RVD-E is engineered to simultaneously elute material from swabs, bind a variety of reverse transcriptase quantitative polymerase chain reaction (RT-qPCR) / PCR inhibitors found in clinical samples, lyse specimens and stabilize nucleic acids in a manner that's compatible with RT-qPCR / PCR. The product consists of a proprietary mixture of peptides, salts, stabilizers, buffers and sodium azide to achieve this task. Reagent RVD-E allows for extraction-free amplification of nucleic acids without performing extractions, centrifugations or other sample manipulations, which may introduce errors, contaminants and/or skew the representation of RNA fragments.

## 02 PRODUCT SIZE

Catalog Number	Volume
7833001	1 mL
7833025	25 mL
7833100	100 mL
7833500	500 mL

## 03 STORAGE & STABILITY

PCR<sub>opsis</sub><sup>™</sup> Reagent RVD-E is shipped and stored at room temperature. The recommended storage temperature is: 4°C ~ 25°C

## 04 OVERVIEW OF PROTOCOL



- 1) Add ~100  $\mu$ L Reagent RVD-E to tube with test swab
- 2) Vortex for ~30 seconds, 3 times, to elute specimen off swab



Reagent RVD-E



Add ~50  $\mu$ L of eluted sample to a thin-walled tube or plate



Heat eluted sample to 95°C as indicated below

Mammalian: 5 minutes  
Viruses: 10 minutes  
Bacteria: 15 minutes

**NOTE:** heating for a longer period of time does not negatively affect results

**NOTE:** samples can be heated in a thermal cycler or heating block

2021 All rights reserved.



Use processed sample into your desired RT-qPCR / PCR mixture

## 05 WRITTEN PROTOCOL

1. Thoroughly mix Reagent RVD-E to ensure homogeneity, but avoid creating bubbles unnecessarily
  1. Reagent RVD-E has a hazy, white color when homogenized and normal settlement occurs if not regularly mixed
2. Elute material from swab:
  1. Add ~100 µL of Reagent RVD-E to transport tube with swab
  2. Make sure the swab is at least partially submerged into Reagent RVD-E
  3. Vortex for 30 seconds, 3 times, to elute sample
3. Specimen lysis & nucleic acid stabilization:
  1. Transfer ~50 µL of eluted sample into a thin-walled PCR tube / plate and then cap tube or apply plate sealer to plate to prevent evaporation
  2. Heat at 95°C
    1. Mammalian: 5 minutes
    2. Viruses: 10 minutes
    3. Bacteria: 15 minutes
    4. NOTE: heating for a longer period of time does not negatively affect results
  3. Let cool at room temperature for ~10 seconds before continuing

## 05 WRITTEN PROTOCOL

4. Pipette up & down to ensure complete mixing
5. Use lysed / stabilized sample in your desired RT-qPCR / qPCR procedure
  1. Lysed / stabilized sample should represent 15%~30% of your final RT-qPCR mixture (i.e., 3~6  $\mu\text{L}$  sample into a total volume of 20  $\mu\text{L}$ ) depending on the polymerase used
  2. You might observe increasing PCR inhibition when your PCR mixture consist of >35% processed sample

## 05 WRITTEN PROTOCOL

Suggested thermocycler parameters for RT-PCR / PCR:

1. Reverse transcription:
  - a. 45°C for 15 minutes
  - b. 95°C for 2 minutes
  
2. PCR amplification (~40 cycles):
  - a. 95°C for 5 seconds
  - b. 55°C for 15 seconds
  - c. 72°C for 15 seconds
  
3. Hold: 4°C

### **NOTE:**

- For most applications, a 3-step PCR amplification set-up is recommended over a 2-step PCR amplification set-up
- The suggested cycles, temperature, and heating times mentioned above may be optimized by the user as needed

## 06 STEP-BY-STEP PROTOCOL WITH FIGURES

### Step 1



Gently invert  
Reagent RVD-E to ensure  
homogeneity

PCRopsis™ Reagent RVD-E

### Step 2



Add ~100  $\mu$ L Reagent RVD-E  
to transport tube containing  
test swab

2021 All rights reserved.

### Step 3



Vortex tube with swab  
and Reagent RVD-E to  
elute specimen off swab

7 of 11

## 06 STEP-BY-STEP PROTOCOL WITH FIGURES


### Step 4



Transfer ~50  $\mu$ L of eluted sample to thin-walled plate or tube

PCRopsis™ Reagent RVD-E

### Step 5



Seal 96-well plate with a plate sealer to prevent evaporation

2021 All rights reserved.

### Step 6

Pre-heated thermal cycler or heating block before applying plate or tube



Recommended heating times:

- Mammalian: 5 minutes
- Viruses: 10 minutes
- Bacteria: 15 minutes

NOTE: heating for a longer period of time does not negatively affect results

Heat eluted mixture at 95°C for recommended time



## 06 STEP-BY-STEP PROTOCOL WITH FIGURES


### Step 7



Add 15  $\mu$ L RT-qPCR mix  
from your desired vendor  
to a new plate

PCRopsis™ Reagent RVD-E


### Step 8



Mix 5  $\mu$ L of heated, eluted  
sample with your  
RT-qPCR mix

2021 All rights reserved.

### Step 9



Detect amplification of  
target genes using your  
desired qPCR equipment

9 of 11

## 07 TROUBLESHOOTING & SUGGESTIONS

1. Reagent RVD-E is optimized for the amplification of gene targets from specimens on swabs and may not be applicable for other applications.
2. For best results, use recently collected samples that have been properly stored since collection.
3. A complete validation study is warranted when using small synthetic swabs and lower volumes of Reagent RVD-E to improve assay sensitivity.
4. Ct cut-offs for assays should be increased, often times to 40 cycles.
5. Ensure that the processed sample consist of <30% of the total PCR mixture, since high concentrations of processed sample may inhibit PCR for some applications.
6. Take care in maintaining the sterility of your Reagent RVD-E stock after use.
7. Heat Reagent RVD-E / sample mixture for a few minutes longer if you observe suboptimal lysis.
8. It's recommended to use the heated Reagent RVD-E + sample mixture for downstream applications within a day.

## 08 CONTACT

Contact our research team if assistance with Reagent RVD-E is necessary ([info@entopsis.com](mailto:info@entopsis.com)). We will try our best to assist with non-intended applications of this product or direct you to alternative products. Any business related questions should be directed to: [Sales@PCRopsis.com](mailto:Sales@PCRopsis.com).



7600 NW 69<sup>th</sup> Ave, Medley, FL 33166

[www.PCRopsis.com](http://www.PCRopsis.com)

[info@entopsis.com](mailto:info@entopsis.com)

[sales@PCRopsis.com](mailto:sales@PCRopsis.com)

Toll-free 1-888-407-5070

PCRopsis™ is a trademark of Entopsis

**NOT FOR RESALE**

v.20210902

PCRopsis™ Reagent RVD-E



2021 All rights reserved.

11 of 11