

# PCRopsis™ Reagent SRVD

## For In Vitro Diagnostic Use

REF #: 7783025, 7783100, 77831000

REF #: 7783001 (for validation purposes only)

Store at room temperature

### **INTENDED USE**

PCRopsis™ Reagent SRVD is intended for extraction-free amplification of RNA or DNA from properly collected and transported saliva samples.

## 01 INTRODUCTION

PCRopsis™ Reagent SRVD is engineered to simultaneously bind a variety of reverse transcriptase quantitative polymerase chain reaction (RT-qPCR) and PCR inhibitors found in saliva specimens, lyse specimens, and stabilize nucleic acids in a manner that's compatible with RT-qPCR / PCR. The product consists of a proprietary mixture of peptides, salts, stabilizers, buffers, RVD Enhancer, and sodium azide to achieve these tasks. Reagent SRVD allows for extraction-free amplification of RNA / DNA from saliva specimens without performing nucleic acid isolation, centrifugations or other sample manipulations, which may introduce errors, contaminants and/or skew the representation of RNA fragments.

## 02 PRODUCT SIZE

| Catalog Number | Volume  |
|----------------|---------|
| 7783001        | 1 mL    |
| 7783025        | 25 mL   |
| 7783100        | 100 mL  |
| 77831000       | 1000 mL |

## 03 STORAGE & STABILITY

PCRopsis™ Reagent SRVD is shipped and stored at room temperature. The recommended storage temperature is: 4°C ~ 25°C

## 04 TRANSPORT MEDIUM COMPATIBILITY

### **RECOMMENDED:**

- no transport medium
- Universal Transport Mediums
- Viral Transport Mediums
- 10 mM Tris pH 8 + antibiotics

### **NOT RECOMMENDED:**

- ◆ Mediums containing guanidinium thiocyanate

NOTE: the user must confirm the compatibility of Reagent SRVD with desired saliva samples

## 05 OVERVIEW OF PROTOCOL



**Reagent  
SRVD**

Mix 20  $\mu$ L Reagent  
SRVD + 20  $\mu$ L  
Sample / well

Heat Reagent SRVD +  
Sample mixture to 95°C for  
recommended time

Use processed sample into  
your desired  
RT-qPCR / PCR mixture  
(sample = <30% of total  
volume)



**Saliva Sample**

NOTE: samples can be  
heated in a thermal cycler  
or heating block

Recommended heating times:

- Mammalian: 5 minutes
- Viruses: 10 minutes
- Bacteria: 15 minutes

NOTE: heating for a longer period of time  
does not negatively affect results

## 06 WRITTEN PROTOCOL

1. Thoroughly mix Reagent SRVD to ensure homogeneity, but avoid creating bubbles unnecessarily
  1. Reagent SRVD has a hazy, white color when homogenized and normal settlement occurs if not regularly mixed
2. Mix 1 volume of saliva sample (20  $\mu$ L) with 1 volume of Reagent SRVD (20  $\mu$ L) in a thin walled tube (0.2 ~ 0.6 mL) or 96-well PCR plate
3. Pipette up & down to ensure complete mixing and then cap tube or apply plate sealer to plate to prevent evaporation
4. Heat diluted sample at 95°C and let cool at room temperature for ~10 seconds before continuing
  1. Recommended heating times at 95°C:
    1. Mammalian: 5 minutes
    2. Viruses: 10 minutes
    3. Bacteria: 15 minutes

NOTE: heating for a longer period of time does not negatively affect results
  2. Make sure the heating device has reached the desired temperature before applying sample

## 06 WRITTEN PROTOCOL

3. Sample heating can be performed using a controlled heating block or thermal cycler
5. Mix heated sample and use lysed / stabilized sample in your desired RT-qPCR / PCR procedure
  1. Lysed sample should represent 15%~30% of your final RT-qPCR mixture (i.e., 3 ~ 6  $\mu$ L sample into a total volume of 20  $\mu$ L) depending on the polymerase used

## 06 WRITTEN PROTOCOL

Suggested thermocycler parameters for RT-PCR / PCR:

1. Reverse transcription:
  - a. 45°C for 15 minutes
  - b. 95°C for 2 minutes
  
2. PCR amplification (~40 cycles):
  - a. 95°C for 5 seconds
  - b. 55°C for 15 seconds
  - c. 72°C for 15 seconds
  
3. Hold: 4°C

### **NOTE:**

- For most applications, a 3-step PCR amplification set-up is recommended over a 2-step PCR amplification set-up
- The suggested cycles, temperature, and heating times mentioned above may be optimized by the user as needed

## 07 STEP-BY-STEP PROTOCOL WITH FIGURES

### Step 1



Gently invert  
Reagent SRVD to ensure  
**homogeneity**

PCRopsis™ Reagent SRVD

### Step 2



Add Reagent SRVD to  
reservoir

### Step 3



Add 20  $\mu$ L of  
Reagent SRVD to wells  
in a 96-well PCR plate



## 07 STEP-BY-STEP PROTOCOL WITH FIGURES


### Step 4



Mix 20  $\mu$ L of sample to wells containing 20  $\mu$ L of Reagent SRVD

PCRopsis™ Reagent SRVD

### Step 5



Seal 96-well plate with a plate sealer to prevent evaporation

### Step 6

Pre-heated thermal cycler or heating block before applying plate



Recommended heating times:

- Mammalian: 5 minutes
- Viruses: 10 minutes
- Bacteria: 15 minutes

NOTE: heating for a longer period of time does not negatively affect results

Heat Reagent SRVD + sample mixture at 95°C for recommended time

## 07 STEP-BY-STEP PROTOCOL WITH FIGURES


### Step 7



Add 15  $\mu$ L RT-qPCR mix  
from your desired vendor  
to a new plate

PCRopsis™ Reagent SRVD

### Step 8



Mix 5  $\mu$ L of heated sample +  
Reagent SRVD mixture with  
your RT-qPCR mix

### Step 9



Detect amplification of  
target genes using your  
desired qPCR equipment

## 08 TROUBLESHOOTING & SUGGESTIONS

1. Reagent SRVD is optimized for the amplification of gene targets from saliva specimens and may not be applicable for other applications or mediums.
2. For best results, use recently collected samples in compatible transport mediums that have been stored at ~4°C since collection.
3. Samples are diluted 50% with Reagent SRVD. As such, you should expect to observe slightly higher Ct's compared to nucleic acid extracted samples.
4. Ensure that the processed sample consist of <30% of the total PCR mixture, since high concentrations of the processed sample may inhibit PCR for some applications.
5. Ct cut-offs for assays should be increased, often times to 45 cycles.
6. Take care in maintaining the sterility of your Reagent SRVD stock after use.
7. Heat Reagent SRVD / sample mixture for a few minutes longer if you observe suboptimal lysis.
8. It's recommended to use the heated Reagent SRVD + sample mixture for downstream applications within a day, although samples may be stable for months at 4°C or -20°C.

## 09 CONTACT

Contact our research team if assistance with Reagent SRVD is necessary ([info@entopsis.com](mailto:info@entopsis.com)). We will try our best to assist with non-intended applications of this product or direct you to alternative products. Any business related questions should be directed to: [Sales@PCRopsis.com](mailto:Sales@PCRopsis.com).



7600 NW 69<sup>th</sup> Ave, Medley, FL 33166

[www.PCRopsis.com](http://www.PCRopsis.com)

[info@entopsis.com](mailto:info@entopsis.com)

[sales@PCRopsis.com](mailto:sales@PCRopsis.com)

Toll-free 1-888-407-5070

PCRopsis<sup>TM</sup> is a trademark of Entopsis

**NOT FOR RESALE**

v.20210902

PCRopsis<sup>TM</sup> Reagent SRVD

