

## **PCRopsis™ Reagent RVD with RVD Enhancer + PCRopsis™ Support enables direct amplification of bacteria from urine**

Francis Buan Hong Lim<sup>1</sup>, Abhignyan Nagesetti<sup>1</sup>, Kevin Moreno<sup>1</sup>, Agustin Galecio<sup>1</sup>, Ian Cheong<sup>1,2,3,4</sup>, Obdulio Piloto<sup>1</sup>.

### **Affiliations:**

<sup>1</sup> Entopsis, Inc., USA

<sup>2</sup> Temasek Life Sciences Laboratory, Singapore

<sup>3</sup> Department of Biological Sciences, National University of Singapore, Singapore

<sup>4</sup> Pathnova Laboratories Pte Ltd, Singapore

Correspondence to: [info@entopsis.com](mailto:info@entopsis.com)

### **Abstract:**

Urinary tract infections (UTI) are a recurrent problem in nursing homes and other settings. Such infections can lead to serious complications if not properly diagnosed early. Symptomatic UTIs result in at least 7 million visits to outpatient clinics, 1 million visits to emergency departments, and 100,000 hospitalizations annually.<sup>1</sup> People suspected of UTIs undergo PCR testing to determine the presence of abnormal bacteria in their urine. In this report, we show that PCRopsis™ Reagent RVD with RVD Enhancer, in combination with PCRopsis™ Support, enables extraction-free direct PCR from urine samples containing gram positive or negative bacteria. These reagents were able to process urine samples containing less than 50 cfu / mL of bacteria. PCRopsis™ Reagent RVD with RVD Enhancer + PCRopsis™ Support offers a robust, low-cost, and fast alternative to nucleic acid extraction for UTI PCR testing.

## Results:

| URINE  | Staphylococcus aureus |       | Pseudomonas aeruginosa |       |
|--|-----------------------|-------|------------------------|-------|
|  | cfu / mL              | Ct    | cfu / mL               | Ct    |
| Reagent RVD<br>with RVD<br>Enhancer +<br>Support | 420,000               | 20.17 | 730,000                | 20.52 |
|  | 42,000                | 23.05 | 73,000                 | 23.93 |
|  | 4,200                 | 27.19 | 7,300                  | 26.83 |
|  | 420                   | 30.89 | 730                    | 29.88 |
|  | 42                    | 35.22 | 73                     | 33.18 |
|  | 4.2                   |       | 7.3                    | 38.49 |

**Table 1. Reagent RVD with RVD Enhancer + Support mediates direct PCR amplification of urine samples with less than 50 cfu / mL of *S. aureus* and *P. aeruginosa*.** The use of these reagents as an alternative to DNA extraction for PCR-mediated detection of diluted bacteria in urine resulted in Ct values ranging from 20 ~ 39. Each sample was tested in duplicates and the indicated Ct is the average of two readings.

| UVT  | Staphylococcus aureus |       | Pseudomonas aeruginosa |       |
|--|-----------------------|-------|------------------------|-------|
|  | cfu / mL              | Ct    | cfu / mL               | Ct    |
| Reagent RVD<br>with RVD<br>Enhancer +<br>Support | 420,000               | 20.66 | 730,000                | 19.83 |
|  | 42,000                | 23.55 | 73,000                 | 22.79 |
|  | 4,200                 | 27.03 | 7,300                  | 25.39 |
|  | 420                   | 30.44 | 730                    | 29.29 |
|  | 42                    | 34.19 | 73                     | 34.14 |
|  | 4.2                   |       | 7.3                    | 36.52 |

**Table 2. Reagent RVD with RVD Enhancer + Support mediates direct PCR amplification of bacterial samples in transport medium with less than 50 cfu / mL of *S. aureus* and *P. aeruginosa*.** The use of these reagents as an alternative to DNA extraction for PCR-mediated detection of diluted bacteria in transport medium resulted in Ct values ranging from 19 ~ 37. Each sample was tested in duplicates and the indicated Ct is the average of two readings.

## Key Conclusions:

- Reagent RVD with RVD Enhancer, in combination with PCR*opsis*™ Support, mediates amplification of gram positive and negative bacteria in urine and transport medium at concentrations less than 50 cfu / mL.
- UTIs usually involve infections with over 10,000 cfu / mL in centrifuged urine.<sup>1,2</sup> The ability to quantitatively detect bacteria in urine over a range expanding more than 400,000 cfu / mL to less than 50 cfu / mL, offers laboratory personal confidence and flexibility in the PCR test.

- The addition of PCR*opsis*<sup>™</sup> Support to PCR*opsis*<sup>™</sup> Reagent RVD with RVD Enhancer is critical for the amplification of low levels of bacteria in urine. Previous PCR studies with bacteria spiked into urine and PCR*opsis*<sup>™</sup> Reagent RVD with RVD Enhancer alone used as an alternative to nucleic acid extraction showed a limit of detection >1000 cfu / mL (data not shown).

## Methods:

### Materials:

- Microorganisms (ATCC): Staphylococcus aureus, Pseudomonas aeruginosa
- BD<sup>™</sup> Universal Viral Transport (UVT)
- IDT: DNA primers & probes (see Table 1)
- Entopsis:
  - PCR*opsis*<sup>™</sup> Reagent RVD with RVD Enhancer
  - PCR*opsis*<sup>™</sup> Support
  - PCR*opsis*<sup>™</sup> 5x PCR Master Mix
- Stellar Scientific: Thin-walled PCR tubes
- Chai: Open qPCR Thermocycler

**Table 1. Primers and probes used for the detection of indicated microorganisms in universal transport medium**

| Microorganism          | Primer | Probe  | Sequence (5'-3')         |
|------------------------|--------|--------|--------------------------|
| Staphylococcus aureus  | S868F  |        | CCACATGCCTCTAATAATG      |
|                        | S1064R |        | GCGATTTTATTTTCTTTTGTAC   |
|                        |        | S1024P | ATGCCATGCCTCCAATATCGC    |
| Pseudomonas aeruginosa | Pa23F  |        | TCCAAGTTTAAGGTGGTAGGCTG  |
|                        | Pa23R  |        | CTTTTCTTGGAAGCATGGCATC   |
|                        |        | Pa23P  | AGGTAAATCCGGGGTTTCAAGGCC |

### Prepare diluted bacterial samples in urine and UVT:

1. Perform 1:10 serial dilutions of bacteria in urine and UVT
2. Mix thoroughly
3. Plate bacterial dilutions to determine colony-forming units (cfu)

### Studies with PCR*opsis*<sup>™</sup> reagents:

1. Add 2 µL of PCR*opsis*<sup>™</sup> Support to 1 mL PCR*opsis*<sup>™</sup> Reagent RVD with RVD Enhancer, and mix thoroughly
  - a. This mixture is stable for 24 hours

2. Thoroughly mix 20  $\mu\text{L}$  of freshly prepared reagent with 20  $\mu\text{L}$  of urine or UVT samples in a thin walled tube (0.2 mL) and cap tube.
1. Incubate the reagent mixture with bacterial sample for 10 minutes at 95°C and let cool at room temperature for ~10 seconds before continuing
3. Add 5  $\mu\text{L}$  of processed sample to 15  $\mu\text{L}$  of PCR mixture.
  - a. PCR mixture:
    - i. PCR<sub>opsis</sub><sup>™</sup> PCR Master Mix, 5X: 4  $\mu\text{L}$ ,
    - ii. IDT primer / probe: 1.5  $\mu\text{L}$ ,
    - iii. Nuclease-Free Water: 9.5  $\mu\text{L}$ ,
    - iv. Processed sample: 5  $\mu\text{L}$ .
4. Run samples on PCR Thermocycler for 45 cycles.
  - a. DNA Amplification:
    - i. 95°C 5 minutes (initial denaturation)
      1. 95°C 30 seconds
      2. 55°C 30 seconds
      3. 72°C 30 seconds
    - ii. 72°C 60 seconds (final extension)
    - iii. 4°C hold

#### References:

1. Schappert, SM. Ambulatory care visits to physician offices, hospital outpatient departments, and emergency departments: United States, 1997. *Vital Health Stat 13*, 1999, vol. 143 (pg. i-i) 1–39
2. Rowe, TA & Juthani-Mehta, M. Urinary tract infection in older adults. *Aging Health*, 2013, 9(5):10.2217.